Use of High-pressure Gas

– Liquid nitrogen-related matters –

CE (Cold Evaporator)

CE safely stores ultralow-temperature liquefied gas inside the tank and either gasifies it with a liquid-to-gas vaporizer (supplied to the **centralized gas piping system**) or supplies it as a liquid (**liquid collection**).

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

Red
mark

Liquid level 600 mm = 352 m³ remaining = 483 L

Red mark = Alert level: Liquid nitrogen collection is prohibited!

An alert is issued when the liquid level reaches 300 mm = 144 m³ = 198 L.

Both supply to gas pipeline and natural evaporation should also be taken into consideration.

Collection from CE

Valve M9


Building No. 8


Valve B

Building No. 9


Collection from CE

Building No. 8

Building No. 9

Safety First!

1. Wear proper clothing and follow appropriate procedures.

- Wear **protection gloves**.

- Learn from **an experienced person**.

2. Low temperature is dangerous!

- Be careful in handling **the cooled container** as well as liquid nitrogen!

3. Ensure good ventilation during use.

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Form

Liquid Nitrogen User Card

Date	Faculty	Division
Laboratory	Name	Amount used
		L

Sample entry

Liquid Nitrogen User Card

Date	Faculty	Division
09/11/12	Urban Environmental Sciences	Applied Chemistry
Laboratory	Name	Amount used
XX Lab	XXX XXXX	20 L

Write down container size.

*Write down container volume.

“15 L because I poured about 3/4” is not acceptable.¹²

Monthly Totaling of User Card

Original liquid nitrogen user record cards for Feb. 2010 Sent to Prof. Munakata

Faculty	Division	Research field	Laboratory
Urban Environmental Sciences	Applied Chemistry	Energy Device	Kanemura Lab.

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- Write clearly.
- Do not write a previous course name.

Any of the above may lead to miscalculation or incorrect billing!

No “Shoplifting”!!

It is against the rule to fill up your container without filling out the form!

Building Nos. 8 and 9

Procedures for liquid nitrogen collection from CE tank

1. Use the key distributed at the time of user registration to unlock the CE yard.
2. Collect liquid nitrogen.
3. Fill out necessary items in the liquid nitrogen user card.
4. Lock the CE yard.

Things You Should Know Prior to Use (regarding safety and cost)

1. Usage fee
2. Factors contributing to unit price increase
3. Unavoidable loss
4. How to prevent explosion
5. Quickly and cheaply
6. Small containers difficult to handle

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44

[illegible]

2. Factors Contributing to Unit Price Increase

Loss at the time of collection decreases usage efficiency and increases unit price.

Actual liquid N₂ used in the first half of FY2017 (Apr. to Sept.)

	Building No.8,RI	Building No.9,10,13
Amount decreased in CE tank	42,030.0 m ³	28,980.0m ³
Total amount used	23,426.9 m ³	20,602.0m ³
Amount of liquid N ₂ collected	13,630.7 m ³	4,516.8 m ³
Amount of gas used	9796.2 m ³	16,085.2 m ³
Usage efficiency	<u>55.7%</u>	<u>71.1%</u>

(First half of FY2016

56.0%

64.6%)

3. Unavoidable Loss

Liquid-level gauge

Pressure gauge

Pressure: 0.7 MPa = Approximately 7 atm

Liquid nitrogen temperature: 77.35 K @ 1 atm
: 100 K @ 7 atm

To collect liquid N₂ at 77 K, cooling of liquid N₂ by gasification is required.


Note that pressurization of CE to supply pipeline gas increases loss at the time of liquid collection.

4. How to Prevent Explosion

Used liquid nitrogen can be recovered

⇒ Eliminate wasting


⇒ Do not expose contents to atmosphere

5. Quickly and Cheaply...


Reduce loss associated with cooling of the container

Leaving a small amount of liquid N_2 at the container base saves time and money because the container is already cool when you make a second addition


Temperature stability of Locator

The superior insulating property of Locator makes it possible to maintain the temperature of the sample at -190°C or less even when the level of liquid nitrogen in the container is less than 5 cm.

6. Small Containers Difficult to Handle

If N_2 is introduced rapidly from the transfer tube, N_2 does not collect inside the container but instead spews out, so no liquid N_2 is collected however long you wait.
→ For small containers, collection may be facilitated by opening the valve to only approx. 20 to 30° and let N_2 flow out gently.

⇔ For ≥ 10 -L containers, **open the valve all the way!**

(If flow rate is too low, N_2 evaporates and never collects.)

- **transporting liquefied gas
by the elevator**

● Rules for transporting liquefied gas by the elevator

1. Only service elevators should be used (*).
2. As a rule, transport should be done by two or more persons.
3. One person should load the "Do Not Enter" sign (Attached Figure) and the liquefied gas container into the elevator at the departing floor. Another person should wait at the destination floor to unload them.

Elevators that can be used for liquefied gas transportation

Elevators shown in red can be used for liquefied gas transportation
(Elevator shown in yellow can be used only for transportation to the first floor of Bldg. No. 9)